I-BOND NF

- SI

Navodila za uporabo

GB

Directions for use

D

Gebruiksaanweisung

ES

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HR

Uputa za uporabu

CMK

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GB

Cobalt-Chrome based dental casting alloy for metal-ceramic dental restorations, contains no beryllium, nickel and gallium.

I-BOND NF is a Cobalt-Chrome based metal alloy for metal-ceramic dental restorations. **I-BOND NF** contains no nickel, beryllium and gallium and fulfils the recommendations of the standards EN ISO 22674 and EN ISO 9693-1 for dental alloys. The relatively low Vicker´s hardness of 285 (HV 10) leads to good and easy milling and grinding behaviour and gives perfect polished surfaces.

COMPOSITION (mass %)	TECHNICAL DATA (Orientation values)	
Co	63	Type
Cr	24	Density
W	8	Vickers hardness HV 10
Mo	3	Coefficient of thermal expansion
		25 - 500°C
		20 - 600°C
Si	1	Melting interval
Nb	<1	0.2%- Elongation limit
		E-modul
		Ductile yield A5
		Casting temperature

Recommendation for use

Wax-up:

Do the wax-up as usual. The wall thickness in wax must not be thinner than **0.35 mm**. Lead wax sprues indirectly. For wax sprues use round wax wires with **Ø 2.5 mm** for single crowns and **Ø 2.5-3.0 mm** for bridges. For greater prosthetic constructions with more than 4 pontics use a distribution bar with **Ø 4.0 - 4.5 mm**.

Melting and Casting:

For investing use phosphate bonded investments for crown and bridge work **Investest K+B Speed** or **Interfine K+B Speed**. Preheat the investment to about **850 to 900 °C**. Hold the end temperature of investment for a minimum **30 minutes**. Refer to manufacturer´s instructions for use for the casting machines.

For **I-BOND NF** use an individual ceramic crucible to prevent contamination with other alloys. Clean crucible after every use.

When melting by induction heating start casting as soon as the ingots have collapsed and oxide net cracks. For melting by flame rotate the reductive zone of the flame around ingots. Start casting as soon as the bath begins to

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vibrate. Allow the cylinder to cool down slowly to the room temperature and then deflask.

Ceramic:

Use ceramic, suitable for porcelain to metal restoration.

- Sand-blast the surface by use of a pencil-blaster. Blast all of the surface on which the porcelain is to be applied with **Al₂O₃ - Interalex 250 µm** and clean the surface (ultrasonically cleaning, under vapour or boiling water, diethyl acetate).
- If **oxide firing** is carried out to check the surface, it need to be fired on **960 °C** at vacuum for **5-10 minutes**. Frame has to show an even dark grey colour. Oxide layer should be removed with **Al₂O₃ - Interalex 250 µm** and cleaned again (vapour, boiling water, ethyl acetate).
- The opaque is then applied on the surface by a first thin wash firing followed by a second covering opaque.
- Firing has to be done according to the ceramic manufacturer´s instructions.
- After each firing step (dentine bake, build-up and glazing) carry out slow cooling down to 700°C.**
- Remove ceramic mechanically. Hydrofluoric acid (HF) corrodes metal frame.

Finishing:

After firing of the metal-ceramic prosthetic work sand-blast visible surface of the frame with **Al₂O₃ - Interalex 50 micr.** and polish it with suitable grinding and polishing instruments. At the end use **Universal polishing paste for Co-Cr-Mo alloys REF 0460** and polish up to high gloss.

Soldering and Welding:

Soldering before firing of the frame can be carried out with suitable base-metal-solder and high temperature flux. For welding with laser use suitable base-metal welding wires.

Final Safety Notes!

Metal dusts are harmful. Use dust extractor and wear gloves. Consider allergic hypersensitivities to contents of the alloy.

Reciprocal Actions!

In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

Warranty!

Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our product are subject to further development. Therefore alterations in construction and composition are reserved.

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