I-MG

Navodila za uporabo / Directions for use
Gebrauchsanweisung / Instrucciones de uso
Uputa za uporabu / Препорака за употреба
Инструкции по использованию / Návod k použití

CZ
INTERDENT s.r.o.
Foerstrova 12, Strašnice
CZ-10000 Praga
tel. +420/274 814 404
fax +420/274 820 130

HR
INTERDENT d.o.o.
Vinogradski odvojak 2d
HR-10431 Sveta Nedelja
tel. +385 1 3873 644
fax: +385 1 3873 617
mail: interdent@interdent.hr

SRB
INTERDENT d.o.o.
Nehruova 44
SRB-11070 Novi Beograd
tel. +381/11 217 53 74
fax +381/11 3 01 75 33

MK
INTERDENT d.o.o.
Bojmija 1-12M
MK-1000 Skopje
tel. +389/224 628 55
fax +389/224 634 12

SK
INTERDENT SK, s.r.o.
Lišcie údolie 57
SK-84231 Bratislava
tel.č. 00421 2 5440554,
00421 2 5440555
e-mail: interdent@interdent.sk

INTERDENT®

Proizvajalec/Producer/Hersteller
SI - 3000 CELJE - Opekarniška cesta 26
Tel.: +386 (0)3 425-62-00 · Fax: +386 (0)3 425-62-10
E-mail:info@interdent.cc · http://www.interdent.cc

CE 0197
Verzija: 08/2014

SI

Zlitina za ulite baze, brez berilija in niklja.

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GB

Non precious alloy for partial denture base without beryllium and nickel.

I-MG is biocompatible Co- Cr- Mo (without beryllium and nickel) partial denture alloy. It has outstanding resistance to corrosion and high strength. **I-MG** can be polished and laser welded. It meets the correspondence to EN ISO 22674.

COMPOSITION (mass %)		TECHNICAL DATA (Orientation values)	
Co	62.5	Type	5
Cr	29.5	Density	8,2 g/cm³
Mo	5.5	Vickers hardness HV 10	365 (HV 10)
Si	1.4	Melting interval	1295-1345°C
Mn, C, N	<1.0	0.2% - Elongation limit	640 MPa (N /mm²)
		Tensile strength	960 MPa (N /mm²)
		E-modul	ca. 220.000 MPa (N /mm²)
		Ductile yield A5	7,5%
		Casting temperature	1450°C

Recommendation for use

Wax up

Sprues should be with round profile with **Ø 3.5-4mm**, avoid rectangular placing of sprues and direct slap of alloy.

Investing:

Use phosphate-bonded partial-denture investment material **Modelcast or Modelcast S**. Preheating temperature of investment is **950°- 1000°C** and depend on moulding and construction of sprues. At total plates is the end temperature **1050°C**. Holding time of end temperature is **45 - 60 minutes**, and depends on the size of the casting flask and filler of the furnace. We ask you to follow recommendation for use of casting machine manufacturer for casting.

Casting:

For **I-MG** use an individual ceramic crucible to prevent contamination with other alloys. Clean crucible after every use.

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When melting by induction heating start casting as soon as the ingots have collapsed and oxide net cracks. For melting by flame heat rotate the reductive zone of the flame around ingots. Start casting as soon as the bath begins to vibrate. Do not overheat the alloy. After cooling down to ambient temperature, deflask the cast and sand-blasted it under 4 bar pressure with **Al₂O₃ - Interloxom 250 µm**. Be careful with clasp sides and stress breakers.

Polishing:

For finishing use carbide, ceramically bonded stones or hart metal milling tools. Use polishing liquid for electrolyte polishing in dental laboratory. Clasps and fitting parts should be protected before polishing with special varnish that prevent uncontrolled material removing. After treating with stones, milling tools and measuring the fit, the cast should be polished with rubber polisher, and at the end with **Universal polishing paste for Cr-Co-Mo alloys REF 460** up to high gloss.

Soldering and welding:

Soldering with Co based lot and high temperature flux. Laser welding with suitable base-metal welding wires.

Final Safety Notes!

Metal dusts are harmful. Use dust extractor and wear gloves. Consider allergic hypersensitivities to contents of the alloy.

Reciprocal Actions!

In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

Warranty!

Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our product are subject to further development. Therefore alterations in construction and composition are reserved.

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