# USER MANUAL

## Aluminium 4-wheeled rollator

Product code: YOLA AIR

### 1. Characteristics

4-wheeled rollator is made of aluminium covered in matt paint. It has large 12" diameter inflated rear and front wheels, height-adjustable handles finished with anatomical grips, brakes with parking clamp, seat with height adjustable backrest and a bag. The structure is based on a crossover frame and an oval cross-section of tubes, thanks to which the resistance to axial pressure is significantly increased. The rollator is equipped with cane and bottle holder for convenient everyday use.

**The set includes:** frame with seat, front and rear pneumatic wheels, backrest, bag, cane holder, bottle holder, pump, mounting screws (2 pcs).

#### Elements of the rollator:

- A. Handgrips
- B. Parking brake
- C. Height adjustment backrest
- D. Seat with handle
- E. Bag
- F. Height adjustment
- G. Cane and bottle holder\*
- H. Folding mechanism for rear wheels
- I. Mounting screws for front wheels
- J. Rear and front 12" pumped wheels

*\* the bottle holder is located at the other side of the frame, not visible on reference picture*

### 2. Application

The rollator is designed for patients with reduced mobility who need stable support and help while walking. It is also recommended to the elderly as an aid in keeping balance in order to improve walking comfort and ensuring safety from

falling. The use of the walker has a positive effect on the return to activity of people with locomotion and tilting to erect position problems. It guarantees better stabilization, support and security while walking.

### **3. Assembling MOUNTING**

**Step 1.** Take the rollator and its parts out of the box.

**Step 2.** Rearrange the front wheels in the correct position (as shown in the reference picture) and secure them with the mounting screws. Pull the red folding mechanism (H) towards you and lower the rear wheels until you hear a characteristic “click”.

**Step 3.** Insert the handles into the openings on the frame and secure them with the locking pins.

**Step 4.** To fully unfold the rollator, grab it by the seat and pull outward. Press the seat on both sides until it clicks into place.

**Step 5.** Attach the backrest by sliding it into the special mount located under the handles. Attach the bag to the front of the rollator using straps.

**Step 6.** Attach the cane holder and the bottle holder to the frame.

**Step 7.** Check that the rollator is assembled properly, fully unfolded, and that no parts or screws are missing. If you notice any damage or missing parts, contact the seller.

**Rollator is ready to use.**

### **FOLDING**

To fold the rollator, repeat above steps in the reverse order.

### **BACKREST HEIGHT ADJUSTMENT**

To adjust backrest height, unscrew the cap located under the parking brake and move the backrest to desired height. Tighten the cap afterwards.

**Note!** Make sure both sides of the backrest are on the same height.

### **HANDGRIPS HEIGHT ADJUSTMENT**

To adjust the height of the handgrip:

- Press the brake handle down to activate the parking brake and immobilize the rollator.
- Press the locking pin located on the frame and adjust height of the handle. Do the same on the other side.
- Unlock the parking brake afterwards.

**NOTE!** Both handgrips should be set at the same height.

## **BRAKE FUNCTION**

- To stop the rollator, pull the brake handle up, then hold.
- When not in motion or when using the seat, lock the parking brake with the parking clamp by pulling the brake handle down.
- To release the parking brake, pull the handle upward.

## **4. Technical data**

|                             |            |
|-----------------------------|------------|
| Overall width               | 68 cm      |
| Overall length              | 82 cm      |
| Seat width                  | 46 x 25 cm |
| Seat height from the ground | 58 cm      |
| Handgrip height adjustment  | 85-100 cm  |
| Rollator weight             | 11,50 kg   |
| Maximum capacity            | 136 kg     |
| Maximum capacity of the bag | 5 kg       |
| Frame colour                | dark blue  |

## **5. Notes**

**Timago International Group does not bear responsibility of improper use of the product, failure to comply safety regulations and misuse.**

## **6. Maintenance**

- Clean the frame with damp cloth and dry thoroughly.
- Do not use strong polishing agents due to the damage of protecting coatings.
- Clean the dirt from moving parts of the rollator regularly (hinges, clamping screws, locks).

## **Tires**

It is recommended to inflate the tires max. to 2,5 bar (36 PSI). The pressure should be adjusted according to the weight of the user and changes in air temperature. Tires may lose air over long periods of storage or non-use. Regularly check the tire pressure and their wear condition. When the tire tread pattern depth is reduced to 1 mm, replace the tires with new ones. Follow these instructions to inflate the tires:

### Slight lack of air

- Press the tire evenly with your hands so that the tire and rim fit evenly.
- Inflate the wheel to the desired pressure.

#### b. No air

- Inflate the wheel to 30% of the total volume and press the tire evenly with your hands so that the tire and rim fit evenly.
- Inflate the wheel to the desired pressure value.

### 7. Storage and transport conditions

The product should not be exposed to direct sunlight, water or moisture.

### 8. Warranty information

All products distributed by our company are covered by a warranty, the terms of which are described in the warranty card available on our website. Please contact the retailer from whom you purchased the product. Please, keep in mind that for warranty purposes, the proof of purchase (receipt or invoice) should be retained.

### 9. Labels

|                                                                                     |                                                                                      |                                                                                     |                             |
|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-----------------------------|
|   | Reference number                                                                     |   | Manufacturer                |
|  | LOT number                                                                           |  | Manufacture date            |
|  | Serial number                                                                        |  | Medical device              |
|  | Note                                                                                 |  | Please read the instruction |
|  | The manufacturer checked for compliance with essential medical devices requirements. |                                                                                     |                             |

Our complete line of products is available from our distributors. Learn more at: **[www.timago.com](http://www.timago.com)**. Thank you for choosing Timago!

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