RÅMJÖLKS MÄTARE REFRAKTOMETER

BRUKSANVISNING

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Colostrum Meter Refractometer

About this product and its use

The VetProvide Colostrum Meter Refractometer is a handheld measuring instrument for checking the quality of colostrum from mares, cows, sheep, and goats. Colostrum is the first milk a newborn animal receives and is rich in antibodies that protect against infection. The refractometer allows you to quickly and easily measure colostrum quality directly in the stable using the Brix scale, which measures in percentage. The higher the Brix value, the more antibodies the colostrum contains.

Measure colostrum quality as soon as possible after birth and before the newborn has nursed. In horses it is especially important that the foal receives colostrum within the first 6 to 12 hours after birth. The gut's ability to absorb antibodies decreases rapidly and is completely gone within 24 hours.

What the Brix value means

The table below shows how to interpret your reading:

Brix value	Quality	Action
23% or above	Good	Colostrum is considered to contain sufficient antibody levels. No immediate action required.
18 to 22%	Fair	Supplementation with donor colostrum or colostrum replacer should be considered. Contact your vet.
Below 18%	Poor	Colostrum likely contains insufficient antibodies. Contact your vet immediately.

Always consult your vet if you are unsure how to proceed after a low reading.

What you will need

In addition to the refractometer you will need:

- Distilled water for calibration
- A small clean container for the milk sample

- A soft clean cloth or paper towel for cleaning

Before you begin: Calibrate the refractometer

The refractometer must be correctly calibrated to give reliable readings. Calibration is done using distilled water and should always be carried out at an ambient temperature of 20 degrees Celsius. Leave the refractometer in the room for at least 15 to 30 minutes before calibrating so the instrument can adjust to room temperature.

C1

Clean the prism

Clean the prism carefully with a soft clean cloth. Make sure it is dry and free of any residue.

C2

Apply distilled water

Place 2 to 3 drops of distilled water on the prism and close the cover gently so the water spreads evenly without air bubbles.

C3

Check the zero point

Point the refractometer towards a light source and look through the eyepiece. Check that the boundary between blue and white sits exactly at 0% on the Brix scale.

C4

Adjust if needed

If the boundary is not at 0%, use the small screwdriver provided and adjust the calibration screw until the boundary sits exactly at 0%.

C5

Complete the calibration

Wipe away the water with a soft clean cloth. Repeat the calibration if you are unsure of the result.

Temperature and accuracy

Calibration should always be carried out at 20 degrees Celsius. If the ambient temperature differs by more than 5 degrees from the calibration temperature, recalibration is recommended to ensure accuracy.

If the refractometer is equipped with automatic temperature compensation (ATC), the instrument can be used directly at temperatures between 10 and 30 degrees Celsius without manual adjustment. Please check with your supplier whether your model has ATC.

Step-by-step measurement instructions

1

Take a milk sample

Express a few drops of colostrum. Discard the very first drops and collect a fresh sample in a clean container.

The first drops may be contaminated and can give a misleading reading.

2

Apply the sample to the prism

Place 2 to 3 drops of colostrum on the prism and close the cover gently so the liquid spreads evenly without air bubbles.

If air bubbles appear, wipe the prism clean and apply a fresh sample.

3

Read the measurement

Point the refractometer towards a bright window or light source. Look through the eyepiece and adjust it until the image is sharp. Read the Brix percentage value where the boundary between the blue and white fields meets the scale.

4

Interpret the result

Compare your reading with the table in the “What the Brix value means” section above and act accordingly.

5

Clean the refractometer immediately after use

Wipe the prism carefully with a soft clean cloth after every use.

Never submerge the refractometer in water for cleaning. This can damage the optical surfaces and internal components.

Important care and maintenance tips

- Always clean the prism carefully after every use to prevent milk residue from affecting the next reading.
- Store the refractometer in the provided case in a dry place, protected from dust and corrosive gases.
- Avoid dropping or knocking the instrument. It is a precision optical device with sensitive surfaces.
- Never touch or scratch the optical surfaces directly with your fingers.
- Check the calibration regularly, at least once before each measurement session.

When should you contact a vet?

Contact your vet immediately if any of the following apply:

- The colostrum Brix value is below 18%.
- The foal, calf, lamb, or kid has not nursed within two hours of birth.
- The newborn appears weak, is not getting up, or is not actively seeking the udder.
- The Brix value is between 18 and 22% and you are unsure how to proceed.
- You are unsure whether the newborn has received enough colostrum.

If colostrum quality is suspected to be insufficient, the newborn's IgG levels should be tested by a vet within 12 hours of birth to confirm that immune protection has been successfully transferred.

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