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Warranty

In accordance with IKA® warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.
The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

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
Flachschlüssel 25 F *
Flat key 25 F *
Clé plate 25 F *

* Nur für S 25 KV - 25 F - IL!
* Only for S 25 KV - 25 F - IL!
* Uniquement pour S 25 KV - 25 F - IL!

Fig. 1


Fig. 2



Fig. 3

IKA-Werke GmbH & Co. KG
Janke & Kunkel-Str. 10
D-79219 Staufen
Tel. +49 7633 831-0
Fax +49 7633 831-98
sales@ika.de

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
Safety instructions

- 
 - The dispersion tool may heats up during operation. In particular in the processing of hot medium, the tool must be cooled before it can be removed from the drive.
 - When in operation the dispersion tool can become hot. If the dispersion tool is not inserted into the drive flange correctly or the medium is hot, it can become extremely hot.
- Don't touch rotating parts during operation.
 - No liquid may get into the drive unit of the disperser.
 - The work has to be interrupted immediately if you notice unusual noise and/ or emission of liquid from the upper rinsing drilling. The seals must be controlled, and replaced if necessary.
 - While working with the disperser, the user must wear his personal protective equipment, in particular safety goggles.
 - Follow the health and safety regulations accident prevention regulations.
 - Note the operating instructions of the disperser unit.
 - Abrasion of the dispersion equipment or the rotating accessories can get into the medium you are working on.
 - Attention!** Never run the dispersion tool dry. Without cooling by the medium the bearing will be destroyed.
 - The dispersion tool is not suitable for continuous operation.
 - Dispersing causes the medium to heat up.
 - Likewise the optimal dispersion duration and rotating frequency must be determined by attempts.
 - The function of the dispersion tools depends on the condition of the sharp edges of rotor and stator. Abrasive media can round off these edges fast, whereby the disperse effect diminishes.

Product information

The shaft bearing of the dispersion tool consists of a two ball bearings and a single-thrust slide ring seal (SiC).
Content of package: shaft and dispersion head.
The plastic materials used are approved by FDA (Food and Drug Administration Approval).
Use the dispersion tool especially for inline operation with the flow chamber.

Application instructions

Application: Suspension, emulsifying, disagglomeration (inline, air-free, sterile).

Spare parts list

Item	Designation	Item	Designation	Item	Designation
1	Cylindrical pin	9	Rotor	2002	Friction ring seal, compl.
2	Shaft tube	10	O-ring 4.7 x 1.42	3003	Pressure spring, right
4	Deep-groove ball bearing	12	O-ring 9.25 x 1.78	3004	Thrust collar
6	Seeger circlip ring	13	Bearing	3005	Slide ring
8	Stator	2001	Shaft, compl.	3008	Counter fit mounting, compl.

Spare parts diagram see front page.
For orders of spare parts, please specify the type of dispersion tool and the designation of the spare part.

Maintenance and cleaning

Allowed procedures	Methods of sterilisation	Notice
Chemical procedures	By germ - killing dissolvents alcohol, phenol,formalin.... Disinfectant remainders are to be removed with germ-free water.	For ball bearings not allowed. Note the chemical resistance of ball bearings and o-rings.
Hot air	Germ - killing by hot-air with 160 to 190°C (approx. 30 min)	Life span of the bearings shortens.

- For pre-cleaning the dispersion tool, operate it in in a dissolvent for dissolving the residuaries of the substance. Due to the large flow rate, rotor and stator are cleaned to a large exten. The dispersion tool must be divided and cleaned immediately after working, so that the adhering residuaries of the substance don't cultivate unwanted bacterial cultures.
For cleaning purposes disassemble the dispersion tool as described below (items see front page):
- Unscrew the rotor **(9)** from the shaft **(2001)**, see **Fig. 2**:
Hold the rotor **(9)** with the flat key. Screw off the rotor **(9)** by turning the shaft **(2001)** in counterclockwise using the shaft key.
 - Unscrew the stator **(8)** from the shaft tube **(2)** (left-hand thread), see **Fig.3**:
Hold the stator with the flat key. Screw off the stator by turning the shaft tube **(2)** in clockwise using the shaft key.
In doing so, the dispersion tool and the stator should be held up in vertical position, thus preventing the spring of the slide ring sealing from falling out.
 - Tap the shaft **(2001)** and the ball bearing **(4)** downward in the direction of the stator thread out of the shaft tube.
 - In this position the parts of the slide ring seal can be removed from the shaft **(2001)**.
 - The slide rings and O-rings must be replaced in the case of wear. The sliding surfaces of the slide and counter ring must be clean.

The assembly of the dispersion tool takes place in reverse order. Make certain that the O-ring **(10)** and the thrust collar **(3004)** are properly fitted during assembly.

Attention! The rotor must not be tightened too strongly, because otherwise the thread will be damaged.

Technical Data

		S 25 KV - 25 G - IL	S 25 KV - 25 F - IL
Working range		Inline	
Stator / Rotor (Ø)	mm	25 / 17	25 / 18
Gap between rotor and stator	mm	0.5	
max. allowable speed	rpm	25000	
max. circumferential speed	m/s	22.2	23.6
Min./max. immersion depth	mm	40 / 85	
Material in contact with medium		FFPM / SiC, AISI 316L	
pH range		2–13	
Suitable for solvents		Yes	
Suitable for abrasive substances		No	
Max. working temperature	°C	220	
Sterilization methods		all meththods	
Working range vacuum	mbar	1	
Working range pressure	bar	6	
Ultimate fineness, suspensions	µm	15–50	5–25
Ultimate fineness, emulsions	µm	1–10	1–5

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