IKA

designed for scientists

IKA VACSTAR digital

Operating instructions

Source language: German

EN 5

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ZH 19

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KO 33

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JA 47

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Fig. 1

Ref.	Description	EN
E	Handle securing screw	
F	Vacuum controller VC 10 stand hole	
G	Stand securing screw	
H	Fan / ventilation slit	
I	Hose connection for suction line d= 8 mm (INLET)	
J	Accessory fastening (emission condenser)	
K	Connection for vacuum controller VC 10 (Mini DIN)	
L	USB interface	
M	RS 232 interface	
N	Water valve connection	
O	Hose connection for pressure line d= 8 mm (OUTLET)	
P	Power supply cable connection	

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Contents

	Page
Device setup	2/3
Declaration of conformity	5
Explication of warning symbols	5
Warranty	5
Safety instructions.....	6
Unpacking	7
Correct use	8
Useful information	8
Setting up.....	9
Commissioning	12
Interfaces and outputs	13
Maintenance and cleaning.....	14
Error codes	15
Accessories	15
Product contact parts.....	15
Technical data.....	16

Declaration of conformity

We declare under our sole responsibility that the product to which this declaration relates is in conformity with directives 2014/35/EU, 2006/42/EC, 2014/30/EU and 2011/65/EU and conforms with the following standards or normative documents: EN 61010-1, EN 60529, EN 61326-1 and EN ISO 12100.
A copy of the complete EU Declaration of Conformity can be requested at sales@ika.com.

Explication of warning symbols

- DANGER

Indicates an (extremely) hazardous situation, which, if not avoided, will result in death, serious injury.
- WARNING

Indicates a potentially hazardous situation, which, if not avoided, can result in death, serious injury.
- CAUTION

Indicates a potentially hazardous situation, which, if not avoided, can result in injury.
- NOTE

Indicates practices which, if not avoided, can result in equipment damage.

Warranty

According to **IKA**'s Terms and Conditions of sale and delivery, this product is covered by a warranty for a period of 24 months. Please contact your dealer for any warranty claims. If you wish, you can send the device directly to our factory. Please provide the sales invoice and state the reasons for your guarantee claim. In this case, you are responsible for shipping and handling costs.

The warranty does not cover wearing parts, nor defects that are the result of improper use, insufficient care and maintenance or failure to follow the instructions in this operating manual.

Safety instructions

For your protection

NOTE

Please read the instruction manual in full before use and follow the safety instructions.

- Keep the instruction manual in a place where it can be accessed easily.
- Ensure that only trained staff use the device.
- Be sure to comply with all safety instructions, directives and all matters of health, safety and accident prevention in the workplace.
- The device and all parts of the device must not be used on people or animals.

DANGER

Always wear personal protective equipment in accordance with the hazard class of the media being worked with - otherwise there are dangers due to spraying liquids, and release of toxic or flammable gases.

- Do not expose human or animal body parts to vacuum.
- Do not work with the device underwater or underground.

Device configuration

DANGER

The vacuum pump IKA VACSTAR digital is not designed to be set up in potentially explosive atmospheres.

- Set up the device in accordance with the chapter “Setting up” and connect the connection lines and interfaces as described.
- Set up the device on a stable, even, non-flammable surface.

- Never work with a faulty or incorrectly connected device.
- Set up the vacuum pumps in appropriate and functioning extractor hoods, or set up an exhaust gas line to them. Ensure that the exhaust gas line cannot be kinked! The maximum permitted exhaust gas line length is 2 metres.
- Prevent occurrence of explosive mixtures; if appropriate connect inert gas for ventilation and/or thinning.

Accessories

- Safe operation can only be ensured when working with accessories as described in the “Accessories” section.
- Carefully observe the operating instructions for additional devices (e.g. rotary evaporators, vacuum controllers), with which the vacuum pump IKA VACSTAR digital is operated.
- The pressure at the gas inlet and outlet must not exceed 1100 mbar.

- Elastic elements may be pressed together under vacuum.
- Only use flexible hose lines.
- Observe your emergency measures for power failure and ensure that the device is put in a safe state (see the chapter: Commissioning, Operating mode).

Work with the device

DANGER

You can prevent the release of solvent vapours into the atmosphere using a downstream emission condenser.

The vacuum pump is not suitable for use with self-igniting materials, materials that are flammable without an air supply, or explosive materials.

Do not operate the pump while it is open. Otherwise there is a danger of severe injuries as a hand may be pulled into the moving parts.

WARNING

Inhaling or coming into contact with media such as poisonous liquids, gases, spray mists, vapours, dusts or biological materials can endanger the health of the user. Ensure that all connections are well sealed and free of leaks if you are working with such media.

- The vacuum pump IKA VACSTAR digital must only be operated under the conditions described in the chapter “Technical data”.
- Prevent release of the materials listed above. Take measures to protect staff and the environment.

- Pay attention to possible interactions or chemical or physical reactions when working with media at reduced pressure and increased temperature.
- There can be electrostatic processes between the medium and the device which can lead to direct danger.
- Some medium may be released due to the residual leakage rate of the device.
- Before commissioning check that all the housing parts are present and fastened to the device.
- Do not lift the pump if the handle is loose and the handle securing screw (E) is loose or missing!
- Only switch the pump on if the pump is standing vertically.
- Connect the hose connections (INLET-OUTLET) and interfaces in accordance with the labelling on the device and the operating instructions.
- Ensure that the temperature of the evacuated medium is below its ignition temperature. The pumping process (compression) increases the temperature of the medium additionally.
- Ensure that vapours containing solvent can be sucked into the pump.
- Do not use the pump to create pressure.
- Observe the permitted pressure at the inlet and outlet side; see the chapter “Technical data”.
- The gas flow must only be regulated/throttled in the suction-side line.
- Use a solenoid valve or a check valve in the suction line if there are several load units.
- When using an emission condenser ensure that the coolant is flowing freely.

Protection of the device

WARNING

Ensure that the substances used with the device are compatible with the materials of the device parts that come into contact with the product; see the chapter “Product contact parts”.

CAUTION

The specified settings on the rating plate must coincide with the actual power supply.

- The device is only disconnected from the power supply when the power or device plug is removed.
- The socket for the power supply cable must be easy to access.
- Covers or parts that can be removed from the device without tools must later be ref

Correct use

Use

⚠ WARNING

The device must not be used for:
- **evacuating biological habitats,**
- **evacuating explosive, corrosive or similar gases,**
- **throughput/use of liquids**

The device is not suitable for:
- **pumping up objects**
- **pumping liquids**
- **use in potentially explosive atmospheres**
- **underwater**
- **underground**

The laboratory diaphragm vacuum pump **IKA VACSTAR** digital is used together with other laboratory devices to create a vacuum. For specific vacuum control, an additional vacuum controller **IKA VC 10** or **IKA** rotary evaporator **RV 10** auto is required.

Standalone operation
The vacuum pump can be operated without a controller. The delivery capacity can be set using the adjustable speed. Generally speaking, the smaller the volume to be evacuated, the lower the speed should be. The final pressure is not controlled.

Operation with accessories
Along with the accessories recommended by **IKA**, the vacuum pump **IKA VACSTAR** digital is suitable for controlled evacuation to a specified final pressure e.g. for operating an **IKA** rotary evaporator. It can also be used for classical separation, filtration or drying in the laboratory.
With an appropriate vacuum controller (e.g. VC 10) or with a rotary evaporator (e.g. RV 10 auto), the VACSTAR digital can be operated automatically, for automatic boiling point detection, scheduled pressure-time curves or programs from the solvent library, for example.
Intended use: Tabletop device

Area of use (only indoors)

The device is suitable for use in residential areas and all other areas.
The safety of the user cannot be ensured:
• if the device is used in conjunctions with accessories not made or recommended by the manufacturer!

- if modifications are made to the device or parts of the device by third parties!
- if the device is operated in contravention of the safety instructions!

Useful information

⚠ DANGER

You can prevent the release of solvent vapours into the atmosphere using a **downstream emission condenser**.

⚠ CAUTION

Special care must be taken when working with gas mixtures that include condensable gases (e.g. solvents). If diaphragm pumps are used with such gases they condense in the pump's air flow. This leads to a build-up of pressure in the vacuum chamber, which in turn influences the effectiveness and service life of the diaphragms and valves.

If liquid or condensate enters the vacuum chambers, the specified minimum working pressure is no longer achieved.

To protect the inner valves and diaphragms against liquid, use an upstream condensate separator, such as a Woulff bottle.

⚠ CAUTION

It is generally recommended that you run the pump with no load at the end of the process (no-load operation, approx. 3 up to max. 5 minutes) to dispel any residues in the pump.

If applicable, vent the recipient. To do this, remove the suction line during manual operation.

⚠ NOTE

Observe the operating instructions of the relevant devices.

Diaphragm vacuum pumps are oscillating displacement pumps. The diaphragms expand the volume of the vacuum chamber to produce an underpressure, which sucks air into the chamber. The air is then pressed out of the vacuum chamber by subsequent compression. Changes between the inlet and outlet are controlled by diaphragm valves.
The physical unit for pressure is 1 Pascal [Pa]. However people often still use the unit millibar [mbar] colloquially.
1 mbar = 100 Pa
1 bar = 10⁵ Pa
1 Pa = 10⁻⁵ bar
1 Pa = 0.01 mbar

A distinction is made between different ranges in vacuum technology.

Rough vacuum:	10 ⁵ to 10 ² Pa	(1000 to 1 mbar)
Fine vacuum:	10 ² to 10 ⁻¹ Pa	(1 to 10 ⁻³ mbar)
High vacuum:	10 ⁻¹ to 10 ⁻⁵ Pa	(10 ⁻³ to 10 ⁻⁷ mbar)
Ultra high vacuum	< 10 ⁻⁵ Pa	(< 10 ⁻⁷ mbar)

The periodic movements of the diaphragms create a vacuum at the pump inlet. Diaphragm vacuum pumps work with a rough vacuum of up to 10² Pa. Lower pressures cannot be achieved with these pumps.
To create fine, high and ultra high vacuums, oil-sealed vacuum pumps, screw pumps, diffusion pumps or turbomolecular pumps must be used

Connection of interfaces

NOTE

Observe the relevant connections (see Fig. 1).

O: Hose connection d= 8 mm OUTLET

Link this connection to the emission condenser with a vacuum hose, or fit a sound absorber to the end.

NOTE

Put the end of the hose in the extractor hood!

Check there is a free outlet on the pressurised side!

Do not use a throttle on the pressurised side and do not close the outlet! Connect the exhaust line to this connection.

K: Connection for vacuum controller VC 10 (Mini DIN)

You can connect the vacuum pump and the vacuum controller VC 10 or the rotary evaporator RV10 auto with the analogue connection cable for precise speed-controlled vacuum control. The vacuum controller detects the pump and switches to speed-vacuum control mode. Two-position control is deactivated. L: USB interface

Connect the vacuum pump IKA VACSTAR digital to a PC with a USB A - USB B cable. Any device software updates can be loaded using the IKA FUT software tool.

M: RS 232 interface

You can connect the IKA VACSTAR digital vacuum pump to a PC using an RS 232 interface cable. The pump can be operated in conjunction with other devices with labworldsoft® laboratory device software. For further information see the chapter "Interfaces and outputs".

N: Water valve connection

Connect the optional water throttle valve RV 10.5001 to the diaphragm vacuum pump. The water flow to the emission condenser is controlled using the water valve. The valve is opened as soon as the pump is switched on.

I: Hose connection for suction line d= 8 mm INLET

Connect the suction line to this connection. Link this connection to the recipient (rotary evaporator condenser, laboratory reactor etc.) with a vacuum hose.

P: Power supply cable connection

Check that the voltage information on the rating plate matches your mains supply. Connect the power supply cable to supply electricity.

IKA VACSTAR digital with vacuum controller VC 10


NOTE

Cannot be operated with optional water valve RV 10.5001 in conjunction with VC 2.4.

Refer to the operating instructions for the vacuum controller VC 10!

When the target value is reached, the suction line is closed by the vacuum valve integrated in the VC 10. The pump speed can be changed by manually setting the speed.

- 1 Recipient (load, e.g. rotary evaporator, reactor)
- 2 Vacuum separator (e.g. Woulff bottle)
- 3 Vent valve
- 4 Pressure sensor
- 5 Vacuum valve/ball valve
- 6 Pump cut-out unit VC 2.4
- 7 Emission condenser

Speed regulation


IKA VACSTAR digital with vacuum controller VC 10 or rotary evaporator RV 10 auto

NOTE

Observe the operating instructions for the devices.

Automatic setting of the operating mode "Speed regulation" as soon as the analogue connecting cable (ref. 8) connects the pump with the vacuum controller or rotary evaporator.

When the target value is reached, the pump speed is automatically reduced to a minimum of "0 rpm". Once the target value is reached, the pump applies suction according to the leakage rate of the vacuum system. Automatic boiling point recognition is possible.

Commissioning

NOTE

Observe the operating instructions for the devices.

Check if the voltage specified on the type plate matches the available mains voltage. If these conditions are met, the device is ready for operation when it is plugged in to the mains.

The socket used must be earthed (grounded contact).

Otherwise safe operation is not guaranteed or the device may be damaged.
Observe the ambient conditions indicated under “Technical data”.

Switching on

- Switch on the device using the main switch (D). All display segments are shown to enable a visual check.

8.8.8.8.

- Software version number display (e.g. vers. 1.00).

0.03

- Operating mode display (A; B). (Factory setting: Operating mode A)

A

The operating mode can be changed when you switch on the device.

Operating mode A:

In this operating mode, the set target value is not saved when the current run comes to an end or the device is switched off. No automatic restart after mains power is cut off.

Operating mode B:

In this operating mode, the set target value is saved when the current run comes to an end or the device is switched off, and the value can be changed.

Operating mode B is required particularly in combination with the pump cut-out unit VC 2.4.

Automatic restart after mains power is cut off.

Changing operating mode:

- Press the rotating/pressing knob (C) and switch on the device at the main switch (D) at the same time to change between the two operating modes.

A

b

In standalone operation or two-position operation in combination with the vacuum controller VC 10 the target speed can be set by turning the rotating/pressing knob (C). Push to start the pump and push again to stop it. A flashing light on the display (B) indicates when the device is running.

Switching on with several interfaces

If the pump is operated via the analogue connecting cable with the vacuum controller VC 10 or the rotary evaporator RV 10 auto, “VC 10” is shown on the display.

VC 10

If the vacuum controller is connected to a PC, “PC” is shown on the display (B).

PC

Switching off

Switch off the device using the main switch (D).

Unplug the appliance from the mains socket to disconnect it from the mains.

Interfaces and outputs

NOTE

Please comply with the system requirements together with the operating instructions and help section included with the software.

The device can be operated in “Remote” mode via an RS 232 or USB interface using the laboratory software labworldsoft®. The RS 232 interface (M) at the back of the device is fitted with a 9-pole SUB-D jack

Accessories: PC 1.1 cable (device to PC)
Required for connecting the 9-pin socket to a PC.


Fig. 10

USB A - USB B cable
Required for connecting the USB interface (L) to a PC.


Fig. 11

Connecting cable
Required for connecting VC 10 and RV 10 auto with the diaphragm vacuum pump **IKA** VACSTAR digital (speed operation).


Fig. 12

Maintenance and cleaning

The device is maintenance-free. It is only subject to the natural wear and tear of components and their statistical failure rate.

Cleaning

For cleaning disconnect the mains plug!
Use only cleaning agents which have been approved by IKA to clean IKA devices.
These are water (with tenside) and isopropanol.
Wear protective gloves during cleaning the devices.
Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.
Do not allow moisture to get into the device when cleaning.
Before using another than the recommended method for cleaning or decontamination, the user must ascertain with IKA that this method does not destroy the device.

Spare parts order

When ordering spare parts, please give:

- device type.
- serial number, see type plate.
- position number and description of spare part, see www.ika.com.

Repair

Please only send devices in for repair that have been cleaned and are free of materials which might present health hazards.

For repair, please request the “**Safety Declaration (Decontamination Certificate)**” from IKA or use the download printout of it from IKA website at www.ika.com.

If your appliance requires repair, return it in its original packaging. Storage packaging is not sufficient when sending the device - also use appropriate transport packaging.

Error codes

If an error occurs, this is shown in an error on the display (B), e.g. Error 4.
Proceed as follows:

- Switch the device off on the device switch,
- Take corrective measure,
- Switch device back on.

Error	Cause	Reason	Solution
Er 3	Overheat	The pump does not run	Switch off the pump, disconnect it from the mains, let it cool down.
Er 4	Motor blocked/overload	The pump does not run	Switch off the pump, disconnect it from the mains Let the pump cool down. Ventilate the vacuum system to normal pressure.
Er 8	Hall sensor defective	The pump does not run	Service
Er 9	BLP self-test/EEPROM error	The pump does not run	Service
Er 10	Outlet Choked	The pump does not run	User needs to check the outlet
Er 11	Membrane stuck	The pump does not run	Air the pump and run at low speed for self-cleaning
Er 12	Connection problem	The pump does not run	Service: check the motor cable connection

In the event that the fault cannot be eliminated using the measures described or another error code is displayed:

- Please contact the Service Department,
- Send the device and a brief description of the fault.

Accessories

- Vacuum controller **IKA** VC 10
- Vacuum control valve **IKA** VCV 1 and VCV 2
- Check valve **IKA** VC 10.300
- Vacuum hose **IKA** VH.SI

Technical data

	Base unit	Value
Delivery capacity max.	m³/h	1.32
Delivery capacity max.	L/min	22
Final pressure	mbar	2
Suction steps		4
Cylinder		4
Connection diameter suction side	mm	8
Connection diameter pressure side	mm	8
Gas ballast valve		no
Inlet pressure min.	mbar	2
Inlet pressure max.	mbar	1030
Boiling point recognition		no
Solvent library		no
Two-position control		yes
Analogue speed-vacuum control		yes
Speed setting option		rotating/pressing knob
Speed min.	rpm	285
Speed max.	rpm	1200
Display		LED
Vacuum sensor		no
Vent valve		no
Condensate separator		no
Emission condenser		no
Volume at minimum pressure	dB(A)	54
Product contact material		PTFE, FFPM, PPS, NBR
Housing material		Cast alu, coated
Appropriate for clean room		no
Explosion-proof		no
Width	mm	150
Depth	mm	370
Height	mm	375
Weight	kg	11.5
Permissible ambient conditions		5 °C to 31 °C at 80 % relative humidity. 32 °C to 40 °C decreasing linearly to a maximum 50 % relative humidity.
Protection class according to DIN EN 60529		IP 20
RS 232 interface		yes
USB interface		yes
Analogue output		no
Voltage	V	100-240
Frequency	Hz	50/60
Power consumption	W	130
Standby power consumption	W	1.5

	Base unit	Value
Protection class		I
Overvoltage category		II
Degree of contamination		2
Use above max. sea level	m	2000
Firmware update		yes
Note		
Protective earth conductor	The protective earth conductor (PE) is not connected to the metal housing. (Internally encapsulated DC power supply)	

Average values.
Subject to technical changes!

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IKA India Private Limited
IKA England LTD.
IKA Vietnam Company Limited