

STEP 1

Close the jaws, then press the "ZERO" switch to set to zero. See figure 1.



FIGURE 1

STEP 2

The skin you use with the skinfold measurement is the belly (about 2.5 cm/1 inch above your right hipbone).



FIGURE 2

STEP 3

While standing, firmly pinch the suprailiac skinfold between your left thumb and forefinger, as shown in Figure 2. Place the jaws of Digital Body Fat Caliper over the skinfold, while you continuing to hold your skind with your left hand, see Figure 3.



FIGURE 3

STEP 4

Press the jaws until you feel a slight pain. Press the "HOLD" switch to keep the measurement. After reading your measurement, press "HOLD" switch to continue. Repeat three times and use the average as your measurement. Refer to the body interpretation chart to determine your body fat percentage and what it means. See Figure 4.



FIGURE 4

SPECIFICATIONS


1. Range: 0-50mm/0 ~ 2 "
2. Resolution: 0.1mm/0.01 "
3. Accuracy: $\pm 0.2\text{mm}/0.01 "$
4. Working temp.: 0 C~+40°C
5. Storage temp.: -20C ~+70°C
6. Battery: 3 V battery

virtuafit

Digital Body Fat Caliper


Body Fat Interpretation Chart for Men

SKINFOLD MEASUREMENT IN MILLIMETERS



Body Fat Interpretation Chart for Women

SKINFOLD MEASUREMENT IN MILLIMETERS



SPECIFICATIONS

- Obtain your body fat measurement in millimeters using the Digital Body Fat Caliper.
- Find where the column with your millimeter reading intersects with the row with your age range.
- The number at this intersection is your body fat percentage.

NOTE

- For measurements over 36mm: Add 0,25% for every millimeter pinched above 36mm.