

USER MANUAL

Aluminium 4-wheeled rollator

Product code: JOHNNIE



1. Characteristics

A four-wheeled, foldable rollator made of an aluminium frame. It is equipped with brakes with a locking mechanism, a seat with a backrest, a crutch/cane holder and shopping bag. It features 8-inch inflatable wheels and height-adjustable handles with anatomical handgrips. The design is based on a cross-braced frame and oval-section tubes, significantly increasing resistance to axial pressure. The rollator is also equipped with suspensions which greatly enhance user comfort both when navigating uneven surfaces (shock absorbing) and while resting in the seat.

The set includes: frame with seat (1 set), handles with brakes (2 sets), backrest (1 pc.), mounting screws with knob (2 sets), shopping bag (1 pc.) crutch/cane holder (1 pc.)

Elements of the rollator:

- A. Handle
- B. Brake
- C. mounting screws with knob
- D. Brake cables
- E. Backrest
- F. Seat
- G. Suspension
- H. Shopping bag
- I. Crutch/cane holder
- J. Upper adjustment of brake
- K. Lower adjustment of brake



2. Application

The rollator is designed for patients with reduced mobility who need stable support and help while walking. It is also recommended to the elderly as an aid in keeping balance, improving walking comfort and ensuring safety from falling. The use of the rollator affects positively on the return to activity of the people with locomotion and tilting to erect position problems. It guarantees better stabilization, support and security while walking.

3. Assembly and use

ASSEMBLY

- To unfold the rollator frame, press down on the sides of the seat so that all parts of the frame snap into place and you hear “click” of the latches under the seat.
- The rollator frame is connected to the handles via brake cables.
Insert the left handle of the rollator into the left opening of the frame. Route the left brake cable forward under the seat of the rollator (a cable holder is located on the rollator frame) and then lead cable to the outside of the frame.
- Lock the handle at the desired height by inserting the mounting screw from the inside of the frame. Tighten the screw from the outside of the frame using the star-shaped knob.
- Perform a similar installation for the right handle.
NOTE! Both handles must be set on the same level.
- Install the backrest by sliding it from top to bottom into the guides located at the front of the rollator handgrips.
- Install the bag by sliding it from top to bottom into the guides located at the front of the rollator frame.
- Install the crutch/cane holder to the lower, horizontal part of the rollator frame so that the top of the holder (the Velcro strap under the handle) is aligned with the bottom part of holder.

UNFOLDING AND FOLDING THE ROLLATOR

- To fold the assembled rollator, pull the red handle located in the centre of the seat. Then fold the frame using the handles.
- To unfold the assembled rollator, unfold the frame using the handles. Then press down on the sides of the seat so that all parts of the frame snap into place and you hear “click” of the latches under the seat.

BRAKES FUNCTION

NOTE! When using the rollator for the first time, you may need assistance from another person, regardless of your physical condition.

- To stop the rollator, pull and hold the brake handle upwards.
- During a stop or while using the rollator seat, lock the brake by pushing the brake handle down.
- To release the brake, pull the brake handle upwards.

HEIGHT ADJUSTMENT OF HANDLES

The height of the handles should be adjusted so that the user can stand upright, as far as he or she is able to, with his or her shoulders in rest position and with slightly bent elbows. It is inadvisable to hang on the rollator with all the body weight.

- Unscrew the star-shaped knob and remove the screw from the rollator frame.
- Set the handle to the desired height.
- Lock the handle at the correct height by inserting the mounting screw from the inside of the frame. Tighten the screw from the outside of the frame using the star-shaped knob.

NOTE! Both handles must be set on the same level.

BRAKES ADJUSTMENT

- **Adjustment of the upper brake:** Unscrew the nut of the cable, pull it up to the appropriate length and lock it by tightening the locking screw.
- **Adjustment of the lower brake:** Loosen the mounting screw with the use of an open-end wrench, pull the black brake cable to adjust its tension, and then, tighten the mounting screw.

4. Technical data

Overall width	64 cm
Overall depth	71 cm
Seat dimensions	44 x 24 cm
Seat height from the ground	55 cm
Handles height adjustment	82-92 cm

Weight	9,0 kg
Maximum capacity	150 kg
Max. capacity of the shopping bag	5 kg
Frame colour	black

5. Notes

Timago International Group does not bear responsibility of improper use of the product, failure to comply safety regulations and misuse.

6. Maintenance

Clean the frame with damp cloth and dry thoroughly.

Do not use strong polishing agents due to the damage of protecting coatings.

Clean the dirt from moving parts of the rollator regularly (hinges, clamping screws, locks).

7. Storage and transport conditions

The product should not be exposed to direct sunlight, water or moisture.

8. Warranty information

All products distributed by our company are covered by a warranty, the terms of which are described in the warranty card available on our website. Please contact the retailer from whom you purchased the product. Please, keep in mind that for warranty purposes, the proof of purchase (receipt or invoice) should be retained.

9. Medical device incidents

Any serious incident involving a medical device must be reported to the manufacturer and to the competent authority of the European Union member state where the user or patient resides.

10. Labels

	Reference number		Manufacturer
---	------------------	---	--------------

	LOT number		UDI code
	Serial number		Medical device
	Note		Please read the instruction
	This product complies with the legal requirements of the European Union regarding medical devices (MDR 2017/745)		

Our complete line of products is available from our distributors. Learn more at: **www.timago.com**. Thank you for choosing Timago!