

MAGRAIN

C L E A R

CONTRAINDICATION

Targets (Patients)

Those who are allergic to the materials used in the pellet and the adhesive should avoid using the product.

PHYSICAL CHARACTERISTICS AND MATERIAL COMPONENTS

Physical Characteristics

1. Pellet: 1 mm in diameter
2. Tape: 7 mm in diameter

Material Components

1. Pellet: iron, chromium, manganese, sulfur, phosphorus, silicon and carbon
2. Tape: acrylic substance

INDICATIONS

MAGRAIN is an acupuncture implement designed to be applied to the skin surface instead of inserting needles in the body.

INSTRUCTIONS

1. Remove perspiration and sebum before application. (Apply the product after bath or skin cleansing.)
2. Peel the product off by twisting the backing sheet at the cut and raising one side of the plaster.
3. Peel one plaster at a time and carefully apply to the affected area.
4. The correct points for application differ with individual complaints. Press therapeutic points (*tsubo*), the stiff, or painful areas with fingertips and apply the product to where the tenderness is the most intense.

PRECAUTION IN RELATION TO THE INSTRUCTIONS

To avoid infection, refrain from using the same plasters for the second time.

PRECAUTIONS

I Precautions (The following patients are advised to use the product with caution.)

1. Those who are prone to rashes should test reaction using one or two plasters before the treatment.
2. Diabetics should consult doctors before the treatment to avoid infection.
3. Those who have tympanic membrane perforation should treat the product carefully as the product may come off and enter the middle (inner) ear when applied to the auricle.

II General Precautions

1. If you feel unwell when using the product, stop using it at once.
2. If any symptoms, such as rashes, reddening, and itching, appear on the skin, stop using the product at once.
3. The condition may seem to worsen in the course of treatment. Remove all plasters and place fewer plasters when symptoms are alleviated.
4. To avoid infection, refrain from scratching the skin after taking plasters off.
5. To avoid rashes, change plasters at least every three days. Put new plasters at different locations for the same affected area.
6. Marks may remain on the skin after plasters are removed. The marks disappear usually in a few days, but in some cases, the discoloration of the skin remains according to the individual physical condition and constitution. Avoid putting plasters on visible parts (e.g. face) when exposed to strong ultraviolet radiation.

III Problems and Adverse Events

1. Problems: falling of pellets and adhesion deficiency
2. Adverse Events: eruption, rubefaction, itch, and wound

IV Other Precautions

1. To prevent accidental ingestion, beware of keeping the product away from children and people suffering from dementia.
2. For pregnant women, application of the product to lower legs may lead to miscarriage and premature delivery.

STORAGE CONDITIONS

1. Plasters should be stored in a cool dry place and protected from the direct sunlight and ultraviolet rays.
2. Do not keep plasters in a refrigerator. The tape quality may deteriorate.

SAKAMURA Lab. & Co.

1 Keshiyama, Kamigamo, Kita-ku, Kyoto 603-8041
JAPAN

Phone: +81-(0)75-701-8171