

Phonak BTE starter kit.

The kit contains:

- 10 wipes
- 1 cleaning spray 30ml
- 2 drying capsules
- 1 drying cup
- 1 cleaning cup
- 1 air blower
- 4 cleaning tablets
- 1 microfiber cloth



Cleaning wipe:

The cleaning wipe is a cleaning product for the maintenance of hearing aids and earmolds. It removes earwax responsible for the decrease in sound quality, increases the life of hearing systems and preserves their effectiveness.



Cleaning spray:

The cleaning spray is a cleaning product for the maintenance of hearing aids and earmolds. It removes earwax responsible for the decrease in sound quality, increases the life of hearing aids and preserves their effectiveness.



Drying capsule:

The drying capsule is indicated for dehumidifying hearing aids and earmolds in order to avoid breakdown due to moisture.



Drying cup:

The drying cup is an essential accessory to be used with drying capsules. It ensures a hermetic environment in order to guarantee efficient dehumidification of hearing aids thanks to the action of dehydrating capsules.



Phonak BTE starter kit.

Cleaning cup:

The cleaning cup is an essential accessory to be used with cleaning tablets. It provides an adequate environment to ensure efficient cleaning of the earmolds through the action of the cleaning tablets.



Air blower:

The airball is designed to remove all water and moisture remaining in the tube of the earmold and in air vents.



Cleaning tablet:

The cleaning tablet is a cleaning product used to clean the earmolds. It removes earwax responsible for the decrease in sound quality, increases the life of the hearing system and preserves their effectiveness.



Microfiber cloth:

The microfiber cloth is exceptionally soft and holds its shape well. The cloth can be used often. It is used for its exceptional ability to absorb water, humidity and residue of cleaning solution. It absorbs eight times its weight in water.

