



designed for scientists



SY stand

/// Data Sheet

The EasySyn reactor systems with double-walled reactor vessels from 500 ml to 5000 ml provide high flexibility using a single, specially designed stand system.

- Quick and easy exchange of different reactor vessels
- Simple alignment of the overhead stirrer
- Adjustment of the operating height to the individual body size.



designed for scientists

Technical Data

Working temperature [°C]	-50 - 200
Attainable vacuum [mbar]	25
Dimensions (W x H x D) [mm]	500 x 1050 x 580