

Chard Advanced Dentistry CAD/CAM Ceramic Characterisation Kit

according to Dr. Simon Chard

At a glance

- Characterization of chairside ceramic restorations
- Consideration of the digital workflow
- Fulfillment of highest aesthetic demands
- Predictably beautiful results

According to
Dr. Simon Chard
Middlesex, UK



CCK01
Chard Advanced Dentistry CAD/CAM Ceramic Characterisation Kit

More Information

Chard Advanced Dentistry CAD/CAM Ceramic Characterisation Kit - Designed for characterisation of chairside ceramic restorations including Empress, Emax and hybrid restorations. This kit has been designed with the aesthetic digital dentist in mind to give

creative control to surface texture, lustre and anatomy after the restoration is milled. This allows the dentist to provide predictably beautiful results if they are just polishing or if they are staining and glazing.

	Fig.	Shank ¹	Ref.-No.	Size	Pieces/Kit
	5936C	HP	807 104 394 504	200	1
	9736G	HP	805 104 042 525	150	1
	882L	HP	806 104 143 524	018	1
	850	HP	806 104 198 524	018	1
	S805	HP	807 104 010 524	018	1
	9770G	HP	803 104 303 533	260	1
	110	HP	100 104 543 000	190	1
	100	HP	090 104 543 000	220	1

¹ 104=HP

	Fig.	Shank ¹	Ref.-No.	Size	Pieces/Kit
	9771G	HP	803 104 000 534	170	1
	9771M	HP	803 104 000 524	170	1
	9771F	HP	803 104 000 514	170	1
	9771C	HP	803 104 000 504	170	1
	DPO06	HP	803 104 372 512	110	1
	943DC	HP	806 104 361 504	100	1
	9734G	HP	805 104 173 525	040	1
	9769M	HP	803 104 000 523	170	1
	9769F	HP	803 104 000 503	170	1

¹ 104=HP